

Thank you

In this unprecedented moment in history, we thank you for your commitment to the well-being of our community, nation, and world at large. Please accept these masks as a token of our gratitude for your immeasurable dedication.

THIS PACKAGE CONTAINS ITEM STYLE: 3-Layer Mask Pocket Ma	ITEMS usk	MADE WITH LOVE BY:
MATERIALS USED BACK LAYER (BODY CONTACTING LAYER):	MIDDLE LAYER: (IF USED)	FRONT LAYER:
☐ 100% Cotton☐ Cotton/Polyester Blend☐ Other:	Pellon 830 Pellon 30 40 810 Other non-woven, non-fusible poly/poly blend interfacing	☐ 100% Cotton ☐ Cotton/Polyester Blend ☐ Other:
MASK CARE BEFORE DISTRIBUTION: Homemade masks should be washed after every use, or as frequently as possible. Wash the homemade mask on a HIGH temperature (preferably at least 158°F or 60°C) standard wash cycle using standard laundry detergent. Dry on high temperature.		
DISCLAIMER Created for Crisis does not recommend the use of these homemade facemasks when surgical masks or N95 respirators are available. These should only be used in non-surgical settings when the healthcare provider (HCP) will be more than 3 feet from a symptomatic patient OR within 3 feet of symptomatic patient		

• Significant exposure to liquid, bodily, or other hazardous fluids

wearing a surgical facemask. HCPs should NOT use these when there is:

- · A high level of infection risk through inhalation exposure
- · A high-intensity heat sources or flammable gases

WARNINGS:

- These homemade facemasks are not surgical masks or Filtering Facepiece Respirators (FFRs).
- · These face masks have not been tested for fluid resistance, particulate or bacteria filtration efficiency, differential pressure, or flammability.

RECOMMENDATIONS:

We recommend monitoring the CDC website below for the latest updates regarding the respirator shortage: https://www.cdc.gov/coronavirus/2019-ncov/hcp/respirators-strategy/crisis-alternate-strategies.html

See www.createdforcrisis.org/labelingandsafety for more information.